

## B-Cool® Wraps and Kits

**NOTE:** Instructions cover multiple product families, images may vary.

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**Intended for single patient use.**



**WARNING:** This product should be used only under the supervision of a medical professional. If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.



**CAUTION:** Never apply a frozen or heated gel pack directly onto the skin.



**CAUTION:** Supervision is strongly advised when applying this product on children and geriatric patients. Patients with circulatory problems should consult a physician or healthcare professional prior to any treatment.



### APPLICATION INSTRUCTIONS

**NOTE:** Before use, follow directions on B-Cool® Gel Packs for Cold or Heat therapy.

1. Open box, remove wrap and gel packs (if included with the kit).
2. Heat or freeze gel packs according to physicians instructions following directions on B-Cool® Gel Pack packaging.
3. Insert prepared gel packs into wrap pockets.
4. Place wrap over affected area and secure with hook.

**NOTE:** Use above pictures to determine placement of straps. Placement of straps will depend on which affected area you will be using cold therapy on and on some applications you may choose to only use one strap.

5. The elastic straps have hook on both ends so they can be easily adjusted as needed for comfort.

**NOTE:** Please refer to liner sew-in label for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product and warranty information or contact a customer care representative at **800-328-1095**.

**NOTICE:** This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

