

Pendulous Abdomen Support with Insert Pocket




Figure 1



Figure 2



Figure 3



Figure 4

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.



WARNING: If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

APPLICATION INSTRUCTIONS

1. Place support pad as low as possible under the tummy, but not so low as to interfere with legs when walking or sitting.
2. Wrap strap around your body and attach the hook to the support. (Figure 1)
3. Grip belt at both ends and stretch forward and secure hook to the support pad. (Figure 2) Adjust as needed until desired level of fit is reached.

NOTE: Optional upper strap can be stretched over the top of the tummy and attach to loop on opposite side for additional support. (Figure 3)

NOTE: For even greater support, lay down on your back and bend knees. Release ends of belt, stretch, and re-attach to the support pad.

NOTE: Please refer to sew-in label on device for care instructions.

Please visit our website at www.birdcronin.com for additional product information or contact a customer care representative at **800-328-1095**.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

**NOT MADE WITH
natural rubber latex**