


# LEVEL-RITE HEEL LIFT

## APPLICATION INSTRUCTIONS



1. Remove heel lift from bag
2. Peel away layers as indicated to provide the desired height
3. Position the lift in the heel area of the shoe

**NOT MADE WITH  
natural rubber latex**

| Reorder No. | Size   | Male Shoe Size | Female Shoe Size | Unit |
|-------------|--------|----------------|------------------|------|
| 0814 0552   | Small  | N/A            | up to 7 1/2      | Each |
| 0814 0553   | Medium | up to 11       | 8 plus           | Each |
| 0815 0554   | Large  | 11 1/2         | N/A              | Each |

### Washing Instructions

Using mild soap, wipe with a warm damp cloth • Air dry

Single Patient Use

DISTRIBUTED BY

**BIRD & CRONIN<sup>®</sup> INC.**

---

---

Quality Medical Products To Help Improve Your Life

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

www.birdcronin.com • ©2015 Bird & Cronin®, Inc. All Rights Reserved

10000555